Liberty Mini

Small Vaginal Probe

Instructions For Use

Please read carefully

TensCare™

INTENDED PURPOSE

The **Liberty Mini** vaginal probe is a medical device designed to be used in the home and healthcare environment to stimulate and train the pelvic floor muscles, helping to treat symptoms of urinary (stress, urge and mixed) incontinence and provide relief from chronic pelvic pain.

The Liberty Mini probe is intended to be used only with muscle stimulation equipment approved to IEC 60601-1, IEC 60601-1-2, IEC 60601-1-11 and IEC 60601-2-10 standards.

INTENDED TARGET POPULATION

The device is intended for use by women aged 18 years and over, in the home and healthcare environment, who suffer from stress, urge or mixed urinary incontinence and chronic pelvic pain, requiring symptomatic relief and management.

TARGET USER PROFILE

The subject device is suitable for use by all who can control the device and understand the instructions.

INDICATIONS AND CLINICAL BENEFITS

The Liberty Mini Probe is intended for use by women aged 18 and over with symptoms of stress, urge or mixed urinary incontinence and chronic pelvic pain. When used as intended, the Liberty Mini Probe can provide the following benefits:

- Reduced frequency of stress incontinence episodes
- Reduced frequency of urge incontinence episodes
- Reduced frequency of mixed incontinence episodes
- Relief from chronic pelvic pain

Pelvic floor muscle training is recommended for improving sexual intimacy. However, it is not intended to treat **any medical issue** and your first step in deciding appropriate therapy should be to consult your professional medical advisor.

Note: Keep the device away from children under the age of 5.

CONTRAINDICATIONS

Follow the Contraindications and Warnings supplied with your stimulator or biofeedback device.

Do not use with electrical stimulation in patients with a history of carcinoma at the site of stimulation.

Do not use in the following circumstances:

- known nickel allergy or sensitivity, due to the risk of allergic reactions or skin irritation;
- pregnancy, unless used under medical advice;
- active vaginal infections or localised lesions;
- impaired sensation in the pelvic region.
- Implanted pacemaker - May affect correct operation of the device.
- Bladder or vaginal infection, or recent surgical scar - Probe insertion may aggravate infection.

PROBE CAUTIONS

Caution: Read all instructions before use. Failure to follow these directions may result in personal injury or damage to the equipment.

Caution: Caution should be used if you have total/partial prolapse of the uterus/vagina because this may impair your ability to fully insert the device and properly position stimulation contacts.

Caution: This probe is intended for single orifice use only. Do not move between anus and vagina. Cross- infection may occur.

Caution: This probe is intended for single patient use only. Do not share your probe with anyone else. Improper treatment or cross-infection may occur.

Caution: It is important that the probe is cleaned before and after each use. *Ineffective cleaning may lead to irritation or infection.*

Caution: Never insert or remove the probe unless the control unit is powered OFF as insertion or removal when stimulation is active may cause discomfort or tissue irritation.

Caution: If tissue irritation occurs, discontinue treatment immediately. *Ask your healthcare professional for advice before continuing further treatment to prevent injury.*

Caution: Do not use a silicone-based lubricant on the metal plates of the probe as it may decrease the effectiveness of the electrical stimulation.

Caution: Electrode surface area is 3.35cm². At high output settings on some devices, this could result in a current density greater than 2mA RMS/cm². High current density may be associated with localised tissue damage or minor non-thermal burns. Increase intensity carefully and reduce if you experience sharp


pain during use or discharge or bleeding after use.

Caution: The stainless steel in the probe's metal plates contain some Nickel. This could cause a reaction if you have a Nickel allergy. *Alternative gold probe specially made without Nickel are available (see X-VPG).*

Warning: Do not use if the pack is opened or damaged.

INSTRUCTIONS FOR USE

1. Connect the lead from the base of the unit to the lead in the probe. Push the pin ends firmly into the pigtail ends of the probe lead.
2. Before using the probe, you will need to visit the toilet.
3. Lubricate the metal electrode surfaces and probe tip with a water-based lubricant, such as TensCare Go Gel.

Caution: Do not use a silicone-based lubricant on the stimulation contacts as it may decrease the effectiveness of the muscle stimulation.

4. Choose a comfortable position, such as lying down on your bed on your side with your knees raised.

Warning: Ensure the unit is switched OFF before inserting the probe.

5. After the wire is securely connected, insert the probe into your vagina, in the same way as you would a tampon, until only the plastic ring at the base of the probe is visible. The metal parts conduct the electrical pulse and should be in contact with the main part of the muscle at all times. The tissues close to the entrance are more sensitive, so you should avoid stimulating them. Probes with long electrodes (the metal part) that run up and down the length of the attachment should always be inserted with the metal parts facing hip-to-hip.
6. After your training session, switch the control unit off, then remove the probe by holding the positioning end ring and gently pulling outwards.

POSSIBLE ADVERSE REACTIONS / SIDE EFFECTS

No serious or long-term adverse effects associated with the use of this device have been reported.

In very rare cases, users may experience mild and temporary adverse reactions, including muscular pain or cramps, vaginal tenderness, increased

vaginal fluid/discharge, vaginal irritation, light vaginal bleeding or urinary tract infection or tingling sensation in the legs.

If any of these symptoms persist, worsen, or cause concern, discontinue use and consult a healthcare professional.

Do not continue using the device if severe pain, heavy bleeding, signs of infection, or any unexpected symptoms occur.

CLEANING

It is important that the probe is cleaned before and after each use. Clean with either an alcohol-free antibacterial wipe such as TensCare Wipes (see X-WIPES) or by wiping with warm soapy water. Rinse under running water and dry thoroughly and return the unit to the storage pouch.

Warnings:

- For single patient use only. May be re-used by the same patient.
- Wash before first use
- If you have active cancer, consult with your clinician before use due to concerns for stimulation potentially spreading cancerous cells.
- Caution should be used if you have suspected or diagnosed heart problems because stimulation devices can affect heart rhythm in certain circumstances.
- Caution should be used if you have suspected or diagnosed epilepsy as unexpected seizures may impair usage.
- Caution should be used if you have a bleeding disorder as stimulation increases blood flow to the stimulated region.
- Do not immerse the probe in a liquid.
- Keep the device away from lint and dust, as long-term exposure to lint or dust may affect the socket.
- Keep the device away from sunlight, as long-term exposure to sunlight may affect the rubber causing it to become less elastic and crack.
- Keep away from children under 5 years of age. Long cable - risk of strangulation in infants.

DISPOSAL

Clean thoroughly and dispose of according to local regulations.

RESIDUAL RISKS

Despite the design, testing, and safety measures applied to this device, some residual risks may remain when it is used as intended.

These include:

- Mild discomfort, irritation, or vaginal sensitivity during or after use
- Muscle soreness or cramps following stimulation
- Tingling sensations, including in the legs
- Minor vaginal bleeding or spotting
- Risk of infection if the probe is not cleaned properly before or after use or if it is shared between users
- Allergic reactions or hypersensitivity to device materials (e.g. nickel or other contact materials)
- Risk of irritation or injury if the probe is inserted incorrectly or excessive force is used
- Risk of burns or discomfort if used with incompatible devices or incorrect stimulation settings
- Reduced or no therapeutic effect if the device is used incorrectly or not as instructed
- Device malfunction or reduced performance if the device is damaged, improperly stored, or used outside recommended conditions.

These risks can be minimised by carefully following all instructions, warnings, and cleaning procedures provided in this manual.

Contact us

If any of these symptoms persist, worsen, or cause concern, discontinue use and consult a healthcare professional. Do not continue use if severe pain, heavy bleeding, signs of infection, or any unexpected symptoms occur.


To report unexpected operation or adverse events, please contact Tenscare Ltd. at:

Email: customer.service@tenscare.co.uk

Phone: +44 1372 723434

TECHNICAL SPECIFICATIONS

Electrical impedance (requirement)	$0.5 \pm 0.2 \Omega$
Material	Electrode body: Plastic ABS PA757 Electrode: Stainless Steel #304 Lead wire: PVC over copper tinsel conductor
Output plug	Fully shielded: touch proof
Probe surface area	3.35 cm^2
Power supply	N/A (the probe doesn't include batteries and does not work by itself. It needs to be connected with a device in order for it to function.)
Weight	24.2 g
Dimensions	86.55 mm (L) x 43.6 mm (P) x 21.5 mm
Service Life	6 months

Applied Part	Vaginal probe
Applied part classification	 Type BF: Equipment providing a degree of protection against electric shock, with isolated applied part. Indicates that this device has conductive contact with the end user.
IP rating	IP22 This medical device is not water resistant and should be protected from liquids. The first number 2: Protected from touch by fingers and objects greater than 12 millimeters. The second number 2: Protected from water spray less than 15 degrees from vertical.
Environmental specifications	
Operating conditions	Temperature range: +5 to 40°C Humidity: 15% to 90% RH non-condensing Atmospheric pressure: 700hPa to 1060hPa
Storage conditions	Temperature range: -25 °C to + 70 °C Humidity: 15% to 93% RH non-condensing Atmospheric Pressure: 700hPa to 1060hPa

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TensCare Ltd.
Advena Ltd.