# Electrodes

Breast 60x130 mm

REF E-CMR60130

## 4 Reusable breast electrode pads for electrotherapy unit

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## **Contents**

### **ENGLISH**

How to Use .....	2
For Toning: .....	2
For Lactation: .....	2
Advice .....	3
Storage life.....	3

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- Electrical Stimulation can build, strengthen, and firm the pectoral muscles giving bulk and lift.
  - Exercise and massage therapy to enhance the size and shape of your breasts
  - Natural methods to enhance your all-natural body
  - Long term results
- All types of sensory stimulation can help initiate lactation. TENS can give discrete sensory stimulation at times when other techniques are not convenient.

### How to Use

Attach to a TENS or EMS stimulator with the lead with red connector (+) closest to the arm, and the black connector (-) closed to the breast bone.

**i** **Note:** The electrode pads must always be used in pairs (two electrode pads on each channel), so that the signal can flow in a circuit.

1) Before use make sure your skin is clean and dry.

**!** **Do not place electrode pads** on skin, which does not have normal sensation or on broken skin.

2) Push the pin ends of the lead wires firmly into the pigtail ends of the electrode pads.

3) Peel the electrode pads from their protective plastic shield by holding and lifting one corner of the pad and pulling.


**!** **Do not** pull on the pigtail wire of the pad.

**Always check** that the unit is **OFF** before applying or removing electrode pads.

4) Apply the adhesive side of the electrode pads as shown in the picture, with the upper and lower pads oriented along the line of the pectoral muscle.

5) Remove the electrode pads from your skin by holding the pad itself and gently pulling.

6) After use always replace the electrode pads on the plastic liner and replace them in the re-sealable plastic bag.

### For Toning:

Use with an EMS unit such as **perfect BEAUTY**, **perfect EMS** or **Sports TENS**.

Choose a muscle toning programme using between 40 and 60 Hz, 250-350  $\mu$ s.

### For Lactation:

Use with a TENS unit such as **perfect TENS** or **perfect mamaTENS**.

For maximum benefit choose programmes with a Burst mode (2 bursts per second) such as **perfect TENS** programme **C** or **perfect mamaTENS** programme **A1**.

- Place the pair of electrode pads on your breasts on each side of the nipple, 1-2 cm away from the areola. Each breast should be connected using lead wires from the same channel: one electrode pad connected to the red lead connector and one to the black.

- Choose the programme and use it for 15 minutes every two hours. Stimulating the sensory nerves in the breast by a TENS machine will imitate a baby suckling (>100 times per minute), which will cause your body to produce Prolactin, the hormone required for lactation. This stimulation should generate a mild tingling sensation. Prolactin is released from the pituitary in about 4-5 discrete pulses over 15-20 minutes with suckling. After this more release does not occur for a while so stimulating the breast for more than 20 minutes does not cause Prolactin release.
- Once you start getting milk in the breast, do not use TENS to replace your regular sessions, because milk staying in the alveoli can turn the alveoli back off. You can use TENS as a supplement stimulation with pumping later or during the night while you sleep to accelerate the changes.

When do you need to induce lactation?

- Adoption
- Maternal-infant separation
- Emergencies/natural disasters

Lactation can be a big struggle for some mums-to-be, especially for adoptive mothers, or mothers who have used surrogates. Breastfeeding does not just mean feeding the baby to fulfil its hunger and achieve the best possible growth, it always means creating a bond between the mother and the baby. Many women desire to induce lactation for the emotional benefits and relational closeness experienced by both mother and baby.

The process of inducing lactation is not complicated but can be slow, and it is recommended to start stimulation about 6 weeks before welcoming the baby. You need to actively stimulate your breasts for 10-15 minutes each at least every 3-4 hours throughout the day. After a few days or weeks (or months) of stimulation you may begin to see your breasts and nipples swelling. They will feel larger and feel sore and under pressure. These are all normal indications that new acini (the milk-producing glandular tissues) are developing inside your breasts. That is necessary before you will be able to produce milk. Once milk appears and you are not yet nursing a child, use a hospital grade breast pump to empty the breast after you use the TENS unit. Emptying the breasts helps you lactate more.

### **Advice**

- The electrode pads supplied are reusable but for single patient use.
- The electrode pads provided are latex-free.
- To obtain the best conductivity through the electrode pads always ensure that they are in good condition and tacky.
- Replace the electrode pads when they lose their stickiness or dry out.
- In very hot weather the gel on the electrode pads may become soft.

### **Storage life**

Electrode pads should last 12 to 20 applications, depending on skin condition and humidity.

Storage life of an unopened pack of electrode pads is 2 years. This may be affected by very high temperatures or very low humidity.

**Important:** Observe all the Contraindications, Cautions and Warnings included with the device that you are using with these electrodes.

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**TensCare Ltd,**
**Advena Ltd.**