
perfectEMS


Instructions for Use
Please read carefully

TensCare™


Quickstart Guide


Connect the device to the electrodes.


Place the electrodes on the area to be treated.


Press and hold the power button to turn the device on.

Select a programme by pressing the P button.

Start the session by increasing the intensity

Press and hold the power button to turn the device off.

Contents

1. Device Description	3
1.1. Introduction	3
1.2. Intended Purpose	3
1.3. Kit Contents	3
1.4. Treatment Options	3
1.5. Controls and Display	5
2. How to Use	6
2.1. Device setup	6
2.2. EMS electrode placement	7
2.3. TENS electrode placement	9
2.4. Completing a treatment	10
2.5. After a treatment	13
2.6. Treatment guidance summary	14
3. Troubleshooting	15
4. Device Care	16
4.1. Cleaning	16
4.2. Storage	16
4.3. Battery replacement	16
4.4. Electrode and lead wire care	17
4.5. Replacement parts/accessories	17
5. Technical Specification	18

1. Device Description

1.1. Introduction

The Perfect EMS is a clinically-proven, dual-channel electrotherapy stimulator with two treatment modes:

- Electrical Muscle Stimulation (EMS) – provides muscle toning, shaping and rehabilitation
- Transcutaneous Electrical Nerve Stimulation (TENS) – provides pain management and relief

The device delivers small electrical signals to the body via body-worn, adhesive electrodes. These signals interact with the muscles and/or nerves. When used as intended, the Perfect EMS can provide the following clinical benefits:

- Reduction in pain
- Improvement in range of motion
- Reduced pain medication reliance
- Improvement in muscle size and strength

1.2. Intended Purpose


Perfect EMS is a medical device designed to be used in the home healthcare environment to improve muscle tone and bulk, as well as provide symptomatic relief and management of acute or chronic localised pain and muscle pain. It is suitable for use by all who can control the device and understand the instructions.

Do not use the device for any purpose other than this intended use.

 **Warning:** Not suitable for use in children without medical supervision.

1.3. Kit Contents

- 1 x Perfect EMS device
- 4 x 50x50mm electrodes (E-CM5050)
- 2 x Lead wires (L-CPT)
- 1 x Detachable belt clip
- 2 x AA batteries
- 1 x Storage pouch
- 1 x Instructions for Use

1.4. Treatment Options

There are two treatment modes to choose from, EMS and TENS.

EMS

EMS is used for muscle training and rehabilitation. The stimulation causes your muscles to contract. The treatment stimulation is intermittent, with 'work' periods in which the stimulation is active and 'rest' periods with no stimulation. This pattern ensures effective muscle engagement and adequate recovery time between contractions.

EMS provides the best results when used to complement physical exercise. Select the desired treatment from the table above.

Prog.	Type	Treatment aim	Muscle Group	Description
A	Preset	Warm up	All	Prepare the muscles before an active training session by stimulating the surface muscle fibres.
B	Preset	Improve endurance	Large (e.g. shoulders, chest, abdominals, glutes, hamstrings, quadriceps)	Stimulate muscles, targeting mainly slow-twitch fibres.
C	Preset	Improve power		Stimulate muscles, targeting mainly fast-twitch fibres.
D	Preset	Improve tone & bulk		Defining muscles in the selected area working mainly on explosive fibres. Grows muscles size.
E	Preset	Improve endurance	Small (e.g. biceps, triceps, calves)	Stimulate muscles, targeting mainly slow-twitch fibres.
F	Preset	Improve power		Stimulate muscles, targeting mainly fast-twitch fibres.
G	Manual	Manual EMS programme	All	Manual EMS programme for using bespoke treatment settings.
H	Manual		All	

TENS

TENS is a clinically proven and safe treatment option for providing pain relief. TENS works in one of two ways:

Pain Gate

Stimulates sensory nerves to block pain signals from being received by the brain. This is achieved through high-frequency settings ranging from 80 to 120Hz. Pain Gate stim provides pain relief whilst being used.

Endorphin Release

Low frequency TENS (typically 2-10Hz) promotes the release of endorphins, the body's natural pain killer. Endorphin Release builds up over 40 minutes and will last several hours after stimulation is ended.

The TENS programmes have either constant or burst output.


Constant output delivers a continuous stimulation, providing a constant, consistent sensation.

Burst output follows an On/Off pattern, giving a pulsing or bursting sensation. This combines the benefits of high-frequency and low-frequency TENS.

Select the desired treatment from the programmes below.

Prog.	Type	Treatment aim	Output	Description
S	Manual	Pain Relief	Constant	Manual TENS programme for using bespoke treatment settings. Recommended for treating acute pain.
T	Manual	Pain Relief	Burst	Manual TENS programme for using bespoke treatment settings. Recommended for treating chronic pain.

1.5. Controls and Display


2. How to Use


2.1. Device setup

Install the batteries

1. Remove the battery cover on the rear of the device.
2. Insert the batteries, ensuring that the positive (+) and negative (-) ends match the corresponding terminals in the compartment.


3. Replace the battery cover.
4. To use the belt clip, slide the belt clip into the rear channels of the device.

Connect the leads

1. Insert the white end of the lead into the channel socket on the rear of the device.

2. Push the metal pins of the lead wires into the sockets on the electrodes.


Position the Electrodes

1. Before applying the electrodes, make sure your skin is clean and dry
2. Select the electrode position based on the selected treatment mode and area of the body being treated. See guidance in the following pages.

The electrodes must always be used in pairs (two electrodes per channel) to ensure the signal can flow in a circuit. The positive (+) electrode should be connected to the red lead wire, and the negative (-) electrode should be connected to the black lead wire.

3. Remove the electrodes from the plastic liner.
4. Apply the electrodes to the body.

i **Note:** always ensure the device is turned off or the intensity is zero before applying the electrodes.


2.2. EMS electrode placement

Electrode placement for EMS is important for obtaining the best results.


Place two electrodes over the bulk of the muscle, with one electrode over the muscle's motor point. The motor point is the area on the skin that is located closest to the motor nerve's entry into the muscle – about 1/3 of the way down the muscle from the spine. Here, it is easiest to trigger a contraction by electrical stimulation.

Experiment by moving the electrodes across the skin until you locate the point over the muscle that gives the cleanest contraction.


Chest


Abdominals


Abdominals


Triceps


Biceps


Shoulders & Traps


Hamstrings


Quadriceps

Calves

Glutes


2.3. TENS electrode placement

The easiest way is to apply the electrodes around/near the source of the pain.


TENS activates the nerves best if it travels along the nerve rather than across it. So, place one electrode further from the spine than the source of the pain, and one closer.

If the pain is in, or close to, your spine, place one electrode either side of the spine.


Neck pain


Upper back pain


Wrist pain


Shoulder pain


Hip pain


Lower back pain


Knee pain


Sciatica

Ankle and foot pain

Elbow pain

2.4. Completing a treatment

Turning on / off

1. Turn the device on/off by pressing and holding the power button located on the top for 3 to 5 seconds.

- i Note:** The device will automatically turn off if the intensity setting is left at zero for more than 5 minutes.

Programme selection

1. Select the desired programme by pressing the P button. (See section 1.4 Treatment options).

Each time you change the programme the intensity level will revert to zero.

When turning on, the device will remember the last programme used.

Treatment time selection

The default treatment time for each programme is 20 minutes.

1. To change the treatment time, press the T button. The 'min' icon should flash on the display.

- Press CH1▲▼ to change the treatment time from either 10, 20, 30, 45, 60, 90 minutes of continuous.
- Press T again to confirm the treatment time.

EMS treatment times should be selected depending on the individual and their training requirements. As with regular exercise, the longer the training session the greater the muscle fatigue and time for recovery. If unsure, start with a 20-minute treatment session and increase/decrease as necessary. TENS can be used for as long as it is necessary. Continuous treatment is fine, but the electrodes should be

repositioned at least every 12 hours to allow the skin to be exposed to the air.

Low frequency TENS works better when the intensity is high enough to cause small muscle contractions. Best results are achieved with a session duration between 20 and 40 minutes. Longer use may cause muscle ache.

Manual programme setting selection


Programmes G, H, S and T are manual programmes that allow bespoke selection of treatment parameters.

The adjustable parameters and related screen icons can be seen below:

Prog	Frequency	Pulse Width	Work	Rest	Ramp
G	1-80Hz	50-350µs	1-40s	0-40s	0-5s
H	1-80Hz	50-350µs	1-40s	0-40s	0-5s
S	1-120Hz	50-250µs	Constant		
T	1-120Hz	50-250µs	2 Hz Burst		

Parameter	Description	Icon
Frequency	The number of pulses per second.	Hz
Pulse Width	The duration of an individual pulse.	µs
Ramp Time	The time required to increase or decrease from zero to the set stimulation intensity, i.e. from work to rest.	▲
Work Time	The time in which the stimulation is active.	W
Rest Time	The time in which the stimulation is at zero.	R
Treatment timer	The duration of the treatment.	min
Constant / Burst Output	The type of stimulation for TENS.	C/B


- To adjust the settings, press the T button. The Hz symbol will flash on the screen, indicating that the frequency setting is ready to be adjusted.

- Press the T button again to cycle through the various adjustable settings.

The corresponding icons will flash, showing which parameter can be modified.

- Adjust the parameters by pressing CH1▲▼.

- After selecting all the parameters and setting the treatment time, press the T button once more to confirm and save your settings.

Starting a session


With the device set-up and the preferred programme and treatment time selected, start the treatment session by increasing the intensity using CH1▲ or CH2▲ (depending on which channel is connected to your body).

Intensity selection

With the treatment started, select the stimulation intensity by using the CH▲▼ buttons.

CH1▲CH1▼ control the intensity of the left-hand lead/electrodes.

CH2▲CH2▼ control the intensity of the right-hand lead/electrodes.

Perfect EMS has 90 levels of intensity. Each press of the ▲ and ▼ buttons will adjust the intensity by 1mA. The intensity can be increased or decreased quickly by pressing and holding the ▲▼ buttons.

Intensity cannot be increased if the device has not been correctly connected to the body (i.e. a closed circuit is made). If the perfect EMS detects a disconnection in either channel it will automatically return the strength to zero in that channel and flash '00'.

To increase intensity, ensure the electrodes are correctly connected to the body and the device.

For EMS programmes, the aim is to produce powerful, visible muscle contractions. The intensity should be increased to roughly 3x the level at which you can first feel tingling, or to as high as you can stand without causing pain.

Use a low intensity for the first session to help you to understand how the machine works. Intensity should be increased

throughout a session and gradually over time. This will help produce effective, pain-free contractions.

For TENS programmes the stimulation sensation should be strong but comfortable. Further increases during use may be necessary if your body becomes used to the sensation. At low intensity levels, it is common for no sensation to be felt. Increase the intensity slowly to find a level that is suitable for you.

Intensity may need to be increased midway through treatment as your body gets used to the sensation. If the sensation becomes uncomfortable, reduce the intensity.

When the yellow LEDs in the channel output sockets light, this indicates the stimulation intensity is high enough to be harmful if not used correctly. Please use with caution.

Keypad lock

If you do not press any keys for 30 seconds, the keypad will lock. This is to avoid accidental changes in setting. To unlock, press either CHI▼ or CH2▼ keys.

Low battery

When the battery is low, the battery icons will display on the screen. The device will automatically turn off after 2 minutes.

Replace with new AA batteries.


2.5. After a treatment

Disconnecting the device

Once the treatment timer reaches zero, the

treatment session will end, and the device turns off automatically.

1. Remove the electrodes from the skin by holding the pad and gently pulling. Do not pull on the lead wires.
2. Replace the electrodes on the plastic liner.
3. Disconnect the leads from the electrodes and device.
4. Place the electrodes in the resealable plastic bag.

5. Return the device, leads and electrodes to the storage pouch.

Frequency of use

Ideally, use EMS programmes once every two days for each muscle with a day of rest in between. It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating one set one day and the other the next day.

Regular use over a period of 4 – 6 weeks should give good initial results, depending on the body shaping, toning, sports training or rehabilitation and recovery needs.

TENS can be used as frequently as necessary and is approved for daily use.


Usage memory

The Perfect EMS tracks the number and duration of use.

To view the usage memory, press T and CH2▼ buttons for 3-5 seconds. The display will show the number of times used and duration in hours used.

Press the same buttons again to return to normal controls.

To reset memory to zero, hold down the T and CHI▼ buttons together for 3 to 5 seconds.

2.6. Treatment guidance summary

	Treatment time	Intensity	Frequency of use
EMS	20 minutes	3x the level at which you can first feel tingling, or to as high as you can stand without causing pain.	Once every two days for each muscle with a day of rest in between.
TENS	High-frequency TENS – as long as necessary to provide pain relief. Low-frequency TENS – 20-40 minutes.	Intensity level which provides strong but comfortable sensation.	TENS can be used as frequently as necessary and is approved for daily use.

3. Troubleshooting

If the device is not working correctly, use the following troubleshooting guide.

Problem	Possible causes	Solution
Display will not turn on	Flat batteries	Replace batteries
	Batteries inserted incorrectly	Remove plastic wrap. Check the positive and negative ends are positioned correctly
	Damaged springs in battery compartment	Contact supplier
Controls will not work	The keypad is locked.	If the key icon is shown, press either CH1- or CH2- to unlock the device.
Low battery display	Low batteries	Replace batteries
No sensation	Intensity is not strong enough	Increase the intensity. Most users will feel something at a setting below 20 mA.
Output will not increase above zero	There is a connection error leading to an open circuit.	Ensure machine is attached to your body correctly.
	Electrodes are not connected to body.	Try using the second lead wire.
	Leads are faulty/damaged.	Purchase replacement if necessary.
Sudden change in sensation	If you stop the stimulation and continue a few minutes later, the signal will feel quite a lot stronger.	Always return the intensity to zero after disconnecting the leads or electrode pads.

If you are unable to resolve the issue, please contact our customer service team:

Email: customer.service@tenscare.co.uk

Phone: +44 1372 723434

 **Note:** The Perfect EMS has no essential performance characteristics. Device failure will not cause unacceptable clinical risk for the user.

4. Device Care

4.1. Cleaning

- Clean the device and lead wires at least once a week by wiping with a damp cloth and solution of mild soap and water. Wipe dry.
- Do not immerse the device in water.
- Do not use any other cleaning solution than soap and water.

4.2. Storage

- Between uses it is important to store the electrodes in their resealable plastic bag with the bag sealed. This will help extend the life of the electrodes.
- Keep the device and all accessories in the storage pouch between uses.

Storage life

- Storage life of an unopened pack of self-adhesive electrode pads is 4 years. This may be affected by very high temperatures or very low humidity.
- The device has no fixed shelf life.
- Unused batteries have a nominal shelf life of 3 years but will usually last longer than this.

4.3. Battery replacement

- When the batteries are low and the battery icon is on the display, replace the batteries.
- Batteries should last at least 15 hours at

50 mA, 300 μ s, 50 Hz.

- The unit will work with rechargeable batteries, but the display may appear dim.

Battery Warnings

- Remove batteries from the device if it is unlikely to be used for a long period. Some types of batteries may leak corrosive fluid.
- Do not pierce, open, disassemble, or use in a humid and/or corrosive environment.
- Do not expose to temperatures over 60°C(140F).
- Do not put, store or leave near sources of heat, in direct strong sunlight, in a high temperature location, in a pressurized container or in a microwave oven.
- Do not immerse in water or get wet.
- Do not short-circuit.
- Do not connect the device unless the battery cover is in place.
- If battery leakage occurs and comes in contact with the skin or eyes, wash thoroughly with lots of water and immediately seek medical attention.
- Keep batteries out of the reach of children to prevent them from swallowing them by mistake. If swallowed by child, contact doctor immediately.
- Never attempt to recharge an alkaline battery. Risk of explosion.

- Do not mix old, new or different types of batteries as this may lead to battery leakage or low battery indication.
- Disposal: Always dispose of batteries and device responsibly according to local government guidelines. Do not throw batteries onto a fire. Risk of explosion.

4.4. Electrode and lead wire care

- The electrodes are reusable but for single patient use. The electrodes should provide up to 20 uses. This may vary depending on skin condition, humidity and storage between uses.
- To maximise the life of the electrodes, always ensure the electrodes are replaced on to their plastic liner and stored and sealed in the plastic bag. When the electrodes dry, this may negatively impact the performance of the treatment.
- The adhesive is the electrode is a peelable hydrogel (water based). To restore the tackiness of the electrode, add a tiny drop of water and spread it around. Avoid adding too much, as it will make the electrode too soft. If this occurs, place the electrodes sticky side up in a refrigerator for a few hours to regain adhesiveness.
- Replace the electrodes when they lose their stickiness. Poor connection may cause discomfort and skin irritation.

- To ensure safety and effectiveness we recommend using TensCare electrodes only.
- The life of the lead wires depends greatly on use. Always handle the leads with care. We recommend replacing the lead wires regularly (about every 6 months).

4.5. Replacement parts/ accessories


The following replacement parts are compatible with the Perfect EMS and may be ordered from TensCare at www.tenscare.co.uk or +44(0) 1372 723434.

E-CM5050	50x50mm electrodes (pack of 4)
E-CM50100	50x100mm electrodes (pack of 4)
E-CR5050	50x50mm Cool Relief infused electrodes (pack of 4)
E-CR50100	50x100mm Cool Relief infused electrodes (pack of 4)
E-TRBG	50x50mm Blue Gel hypoallergenic electrodes (pack of 4)
E-CMR60130	Breast electrodes (pack of 4)
L-CPT	Lead wire
B-AA	AA batteries (pack of 2)

5. Technical Specification

Waveform	Symmetrical Bi-Phasic
Amplitude (over 1 kOhm load)	90 mA zero to peak +ve in 90 steps +/- 10%
Frequency range	1-120Hz
Pulse width range	5-250µs
Output plug	Fully shielded: touch proof
Channels	Dual channel
Power supply	2 x AA alkaline (two AA batteries)
Weight	75g (without batteries)
Dimensions	120 x 60 x 20 mm
Safety Classification	Internal power source. Designed for continuous use. No special moisture protection.
Service life	>5 years
Software version	PFEMS_M6832_R4
Applied part	Skin surface electrodes
Applied part classification	Type BF
IP rating	IP22
Environmental Specifications	
Operating conditions	Temperature range: 5 to 40°C Humidity: 15 to 90% RH non-condensing Atmospheric pressure: 700 hPa to 1060 hPa
Storage conditions	Temperature range: -25 to +60°C Humidity: 15 to 90% RH non-condensing Atmospheric pressure: 700 hPa to 1060 hPa

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Thank you for using Perfect EMS...

At TensCare, we strive to provide the best possible products and services. We value your feedback and continuously work to enhance our offerings. Understanding how our products benefit you is important to us. If you have any insights or suggestions, please get in touch.

Contact us

-  TensCare Ltd.
9 Blenheim Road,
Epsom, Surrey,
KT19 9BE, U.K.
-  +44 (0) 1372 723434
-  www.tenscare.co.uk
-  customer.service@tenscare.co.uk

Follow us

-  TensCare Ltd
-  @tenscareltd
-  @TensCareLtd
-  TensCare Limited

EU Declaration of Conformity

Tenscare Ltd herewith declare that the product meets the provisions of REGULATION (EU) 2017/745 which apply to it. The medical device has been assigned to class IIa according to Annex IV of the REGULATION (EU) 2017/745. The product concerned has been designed and manufactured under a quality management system according to Annex IX, of REGULATION (EU) 2017/745. This EU declaration of conformity is issued under the sole responsibility of the manufacturer. No "Common Specification" is applicable.

www.tenscare.co.uk

TensCare Ltd.

9 Blenheim Road, Epsom,
Surrey, KT19 9BE, U.K.

EU DoC

Complies with Medical Device
Regulation (EU) 2017/745.

Advena Ltd.

Tower Business Centre,
2nd Flr, Tower Street, Swatar,
BKR 4013, Malta

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