
unifit

Instructions for Use
Please read carefully


TensCare™

Quickstart Guide


1

Connect the device to the electrodes.

2

Place the electrodes on the area to be treated.

3

Press and hold the power button to turn the device on.

4

Select the mode by pressing the M button.

5

Select a programme by pressing the P button.

6

Start the session by increasing the intensity.

Contents

1.	Device Description	3
1.1.	Introduction	3
1.2.	Intended purpose	3
1.3.	Kit contents	3
1.4.	Controls and display	4
2.	How to Use	4
2.1.	Device setup	4
2.2.	Completing a treatment	6
2.3.	After a treatment	7
3.	Treatment options	8
3.1.	TENS treatment	8
3.2.	EMS treatment	11
3.3.	Massage programmes	16
4.	Troubleshooting	16
5.	Device Care	17
5.1.	Cleaning	17
5.2.	Storage	17
5.3.	Charging Unifit	18
5.4.	Electrode and lead wire care	19
5.5.	Replacement parts/accessories	19
6.	Technical Specification	20

1. Device Description

1.1. Introduction

Unifit combines the treatment capabilities of a TENS device, an EMS or NMES device, and a relaxing massage stimulator in one unit.

TENS – Transcutaneous Electrical Nerve Stimulation

EMS – Electrical Muscle Stimulation

1.2. Intended purpose

Unifit is a medical device designed to be used in the home healthcare environment to improve muscle tone and bulk, as well as provide symptomatic relief and management of acute or chronic localised pain and muscle pain. It is suitable for use by all who can control the device and understand the instructions.

TENS stands for Transcutaneous Electrical Nerve Stimulation. TENS can be used to provide symptomatic pain relief including chronic pain, postoperative pain and joint pain.

EMS stands for Electrical Muscular Stimulation. EMS is used for muscle training, toning and bulking, as well as rehabilitation post-injury.


Unifit is designed to be used in the home healthcare environment. Do not use the device for any purpose other than it's intended use.


Warning: Not suitable for use in children without medical supervision.

1.3. Kit contents

- 1 x Unifit device
- 4 x 50x50mm electrodes (E-CM5050)
- 2 x Lead wires (L-CPT)
- 1 x Detachable belt clip
- 1 x Charging adapter (UK/EU) (X-UNICH)
- 1 x Storage pouch
- 1 x Instructions for Use

1.4. Controls and display


2. How to Use

2.1. Device setup

Connect the leads

1. Insert the white end of the lead into the channel socket on the rear of the device.

2. Push the metal pins of the lead wires into the sockets on the electrodes.


Position the electrodes

1. Before applying the electrodes, make sure your skin is clean and dry
2. Select the electrode position based on the selected treatment mode and area of the body being treated. See electrode placement guidance for each treatment in section 3.0

The electrodes must always be used in pairs (two electrodes per channel) to ensure the signal can flow in a circuit. The positive (+) electrode should be connected to the red lead wire, and the negative (-) electrode should be connected to the black lead wire.

3. Remove the electrodes from the plastic liner.

4. Apply the electrodes to the body.

i **Note:** always ensure the device is turned off or the intensity is zero before applying the electrodes.

Turning on / off

1. To turn the device on, press and hold the power button for 1 second.
2. To turn the device off, press and hold the power button for 3 seconds.

Note: The device will automatically turn off if the intensity setting is left at zero for more than 5 minutes.

Mode and programme selection

1. Select the desired mode by pressing the M button. (See section 3.0 Treatment options).
2. With the mode chosen, select the programme by pressing the P button. (See section 3.0 Treatment options).

Each time you change the programme the intensity level will revert to zero.

When turning on, the device will remember the last programme used.

Treatment time selection

The default treatment time for EMS preset programmes is shown on the display. For TENS programmes treatment time can be changed.

1. To change the treatment time, press the settings button.

2. Use CH2 ▲▼ to select the preferred treatment time. You can select between 5 - 90 minutes or continuous (marked as C).

3. Press the settings button again to confirm the treatment time.

Manual programme setting selection

When a programme has manual settings available, 'M' will be displayed.

In these programmes you are able to select and change the following parameters:

Frequency (Hz)

Pulse Width (μ s)

Treatment Time (min)

Work Time (EMS programmes only) (s)

Rest Time (EMS programmes only) (s)

Ramp Time (EMS programmes only) (s)

1. To adjust the settings of a programme, press the settings button.
2. Use the CH1 ▲▼ buttons to cycle between the different settings you wish to change. The selected setting will flash.
3. Use CH2 ▲▼ to increase or decrease these settings.
4. To confirm your manual settings, press the settings button again. Alternatively, the device will automatically confirm your settings if no button is pressed for 5 seconds.

Warning: In the EMS manual programmes a warning triangle will be displayed if the Rest period is less than $\text{Rest Time} = \text{Work Time} * (\text{WorkHz} - 16.66)$

/ 16.66 seconds

To reset programmes settings to default, press the settings button to adjust the settings. Then press and hold the mode and programme buttons simultaneously for 3 seconds.

2.2. Completing a treatment

Starting a session

With the device set-up and the preferred programme and treatment time selected, start the treatment session by increasing the intensity using CH1▲ or CH2▲ (depending on which channel is connected to your body).

Intensity selection

Select the stimulation intensity by using the CH▲▼ buttons. CH1▲ CH1▼ control the intensity of the left-hand lead/electrodes. CH2▲ CH2▼ control the intensity of the right-hand lead/electrodes.

The intensity can be increased or decreased quickly by pressing and holding the ▲▼ buttons.

At low intensity levels, it is common for no

sensation to be felt. Continue to increase the intensity until the sensation is strong but comfortable.

During a session, your body can become used to the sensation. If this is the case, increase the intensity levels.

The intensity controls will not operate unless the Unifit has been properly connected to the body. If the Unifit detects a disconnection on either channel, the intensity will rest to 0mA and the connection error icon will flash. This is a safety feature designed to prevent sudden surges in stimulation.

When the yellow LEDs in the channel output sockets light and a white circle is displayed next to the intensity levels, this means the output is high enough to be harmful if not used correctly. Please use with caution.

See guidance on intensity selection for each treatment mode in section 3.0.

Keypad lock

If no key is pressed for 30 seconds, Unifit will auto-lock.

The lock indicator will flash and then remain on.

This is to avoid accidental changes in setting.

To unlock, press either CH1▼ or CH2▼ keys.

Low battery

The battery status can be seen in the top right corner of the display. When the remaining battery life is below 20%, the battery icon will flash.


Please charge your device.

2.3. After a treatment

Disconnecting the device

Once the treatment timer reaches zero, the treatment session will end, and the device turns off automatically.

1. Remove the electrodes from the skin by holding the pad and gently pulling. Do not pull on the lead wires.
2. Replace the electrodes on the plastic liner.
3. Disconnect the leads from the electrodes and device.
4. Place the electrodes in the resealable plastic bag.


5. Return the device, leads and electrodes to the storage pouch.

Usage memory

Unifit tracks the number and duration of use.

To view the usage memory, press and hold the settings and phase skip buttons for 3 seconds. The display will show the number of times used and duration used in hours and minutes.

Press the same buttons again to return to normal controls.

To reset memory to zero, hold mode and programme buttons together for 3 seconds. This can only be done whilst viewing the usage memory.

3. Treatment options

There are three treatment options to choose from; TENS, EMS, and MASSAGE.

3.1. TENS treatment

How it works

TENS is a clinically proven and safe treatment option for providing pain relief. TENS works in one of two ways:

- **Pain Gate**

Stimulates sensory nerves to block pain signals from being received by the brain. This is achieved through high-frequency settings ranging from 80 to 120Hz. Pain Gate stim provides pain relief whilst being used.

- **Endorphin Release**

Low frequency TENS (typically 2-10Hz) promotes the release of endorphins, the body's natural pain killer. Endorphin Release builds up over 40 minutes and will last several hours after stimulation has finished.

Programme selection

The TENS programmes have either constant or burst output.

Constant output delivers a continuous stimulation, providing a constant, consistent sensation.

Burst output follows an On/Off pattern, giving a pulsing or bursting sensation. This combines the benefits of high-frequency and low-frequency TENS.

Programmes with modulated output types can be used to help prevent you getting used to the stimulation and it becoming less effective.

There is not one programme for a particular condition, and the best choice varies from each person, even if they have the same type of pain.

While we provide recommended uses, we advise trying different programs to find the one that works best for you.


Prog.	Type	Output	Type of Relief	Recommended use
1	Preset	Constant	Pain Gate	For first use of TENS. These programmes have slightly different parameters so choose the one that works best for you.
2	Preset	Constant	Pain Gate	
3	Preset	Constant	Endorphin	For chronic pains such as back pain and arthritis
4	Preset	Burst	Pain Gate + Endorphin	For radiating pain in arms and legs and deep muscular pain (Osteoarthritic pain in the knee, Sciatica central pain). These programmes have slightly different parameters so choose the one that works best for you.
5	Preset	Burst	Pain Gate + Endorphin	
6	Preset	HAN	Pain Gate + Endorphin	For maximum production of the two central opioid peptides associated with pain relief.
7	Preset	HAN	Pain Gate + Endorphin	These programmes have slightly different parameters so choose the one that works best for you.
8	Preset	Frequency modulated	Pain Gate + Endorphin	Settings modulate to prevent accommodation occurring.
9	Preset	Frequency and pulse width modulated	Pain Gate + Endorphin	
10	Preset	DTENS	Pain Gate	Treating nerve rich areas (e.g. facial pain, neck pain).
11	Manual	Constant	Pain Gate or Endorphin	For manual settings control.
12	Manual	Burst	Pain Gate +/ or Endorphin	

Electrode placement


The easiest way is to apply the electrodes around/near the source of the pain.

TENS activates the nerves best if it travels along the nerve rather than across it. So, place one electrode further from the spine than the source of the pain, and one closer. If the pain is in, or close to, your spine you can place one electrode either side of the spine.


Neck pain


Shoulder pain


Wrist pain


Lower back pain


Upper back pain


Period pain


Elbow pain


Hip pain


Sciatica

Knee pain

Ankle and foot pain

Intensity selection

The stimulation sensation should be strong but comfortable. Further increases during use may be necessary if your body becomes used to the sensation. At the start of the session, increase the intensity slowly to find a level that is suitable for you.

Frequency of use

It is not recommended to use TENS programmes for longer than 90 minutes.

The low frequency programmes work best when the strength is high enough to cause small muscle contractions. Best results are achieved with a session duration between 20 and 40 minutes. Longer use may cause muscle ache.

3.2. EMS treatment

How it works

EMS is used for muscle training and rehabilitation. The goal of electrical muscle stimulation is to achieve visible muscle contractions. The electrical impulses

delivered from the device stimulate the nerves to send signals to a specifically targeted muscle, which reacts by contracting, just as it does with normal muscular activity.

Mode and programme selection

Each mode consists of three phases, a warm-up phase, and training phase and a cool-down phase. This allows for effective muscle training whilst reducing the chances of injury or fatigue.

The treatment stimulation is intermittent, with 'work' periods in which the stimulation is active and 'active rest' periods where the stimulation is reduced. This pattern ensures effective muscle engagement and adequate recovery time between contractions.

It may improve progress if you combine voluntary contraction with stimulation. For best results, always try to contract voluntarily during the work phases.

There are four EMS modes to choose from:

EMSI	For small muscles such as those in the hands.
EMSII	For medium sized muscles such as arms and calves
EMSIII	For large muscles such as thighs, buttocks and abdomen
EMSIIII	Advanced muscle mass mode for large muscle mass

With the mode selected, choose from one of the following programmes:

Prog	Application	Description
1	Muscle calming	Relax the muscles as much as possible and to promote the body's natural endorphins to promote pain relief and to improve the blood circulation and provide oxygen into the muscle.
2	Resume training	Promotes slow twitch fibres to build muscle strength to help reduce muscle atrophy ready for resuming training activities.
3	Resume training	Progress from programme 2 as tolerance increases.
4	Resistance training	Improving the capacity to develop very high level of muscle force over a long period of time. Improving the efficacy of oxygen consumption at the muscle level and the capacity to withstand toxin accretion, such as lactic acid. For sports activities that require very high levels of prolonged muscle activity: Rowing, Cycling, Middle distance running.

5	Resistance training	Progress from programme 4 as tolerance increases. Alternative application: Lipolysis Increasing blood circulation and modifying the metabolism of the lipocytes. This helps to stimulate the subcutaneous deposits of fat, reducing or eliminating the 'Orange Peel' effect of the skin surface
6	Maximum muscle contraction	To increase muscle bulk and volume and to improve muscle force. Searching for muscular hypertrophy.
7	Muscle toning	Strengthening the muscles, improving blood circulation and capillary bed density. Ideal for applying to the Thigh, Legs, Bottom and Abdomen.
8	Muscle toning	Similar to programme 7, but adding bulk in preference to endurance
9	Fore output, anaerobic activity	Increasing the muscle capacity to a level of instantaneous maximum muscle force, changing muscle force into explosive action. Used for all activities requiring maximum muscle output in a very short space of time, such as Judo, short distance sprinting, throwing the discus or shot-put.
10	Manual	Select the frequency, the pulse width, the work and active rest time. At first use, the default values are shown. Note: There are two manual programmes as the device will remember your previously selected settings. This allows you to create two custom programmes.
11	Manual	

Electrode placement

Place two electrodes over the bulk of the muscle, with one electrode over the muscle's motor point. The motor point is the area on the skin that is located closest to the motor nerve's entry into the muscle – about 1/3 of the way down the muscle from the spine. Here, it is easiest to trigger a contraction by electrical stimulation.

Experiment by moving the electrodes across the skin until you locate the point over the muscle that gives the cleanest contraction.

Note: The +ve electrode needs to be plugged in to the red lead wire connector, the -ve electrode to the black lead wire connector.


Intensity selection

The aim of EMS is to produce powerful,


visible muscle contractions. The intensity should be increased to roughly 3x the level at which you can first feel tingling, or to as high as you can stand without causing pain.

Use a low intensity for the first session to help you to understand how the machine works. Intensity should be increased throughout a session and gradually over time. This will help produce effective, pain-free contractions.


Chest


Biceps


Abdominals


Abdominals


Triceps


Shoulders & Traps


Quadriceps


Hamstrings

Calves

Glutes

Frequency of use

Use EMS programmes once every two days for each muscle with a day of rest in between. It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating one set one day and the other the next day.

Regular use over a period of 4 - 6 weeks should give good initial results.

3.3. Massage programmes

Massage programmes can be used to relieve muscle stiffness. Experiment with the different programmes to find the one that you find the most relaxing.

Programme 1 alternates between two strong squeezing/kneading feelings.

Programme 2 alternates between slow and fast tapping.

Programme 3 alternates between gentle and strong kneading.

Programme 4 has five steps of gradually increasing strength of kneading, alternating with tapping.

Programme 5 & 7 alternates between slow and medium-strong tapping.

Programme 6 is a gentler form of programme 5.

Programme 8 is a pumping sensation roughly equivalent to your heart rate, that slowly increases in speed over the course of the 20 min treatment. It is especially effective with swollen legs and veins.

Programme 9 is squeezing/kneading in 6 steps of increasing intensity.

Programme 10 generates a tapping sensation, gradually speeding up, then slowing down.

Note: There are no restrictions to the frequency of use of massage programmes.

4. Troubleshooting

If the device is not working correctly, use the following troubleshooting guide.

Problem	Possible causes	Solution
No display	Flat battery	Charge battery. See section 5.3
	Damaged battery	Contact supplier. There is no maintenance or calibration other than replacement of gel pads and battery charging. Do not modify the device as this may reduce safety or effectiveness.
Low battery display	Low batteries	Charge battery
Controls will not work	Keypad is locked	If the lock icon is shown, press either of the - buttons to unlock the device.

No sensation	Intensity is not strong enough	Increase the intensity until you feel a strong but comfortable sensation.
Intensity setting returns to 0mA.	There is a connection error leading to an open-circuit.	Ensure the pads are correctly attached to your body and lead wire inserted to both the pads and device.
	Faulty pads or lead wires.	Try using the other lead wire or electrode pads. If problem occurs, it may be necessary to purchase replacements.
Warning icon is flashing	EMS manual mode WORK/REST settings have been adjusted outside the recommended ratio.	Continue treatment with caution and stop if stimulation becomes uncomfortable. Manual programme settings can be reset by pressing and holding the mode and programme buttons simultaneously.

If you are unable to resolve the issue, please contact our customer service team:

Email: customer.service@tenscare.co.uk

Phone: +44 1372 723434

Note: Unifit has no essential performance characteristics. Device failure will not cause unacceptable clinical risk for the user.

5. Device Care

5.1. Cleaning

Clean the device and lead wires at least once a week by wiping with a damp cloth and solution of mild soap and water. Wipe dry.

- Do not immerse the device in water.
- Do not use any other cleaning solution than soap and water.

5.2. Storage

Between uses it is important to store the electrodes in their resealable plastic bag with the bag sealed. This will help extend the life of the electrodes.

Keep the device and all accessories in the storage pouch between uses.

Storage life

- Storage life of an unopened pack of self-adhesive electrode pads is 4 years. This may be affected by very high temperatures or very low humidity.
- The device has no fixed shelf life.
- Unused batteries have a nominal shelf life of 3 years but will usually last longer than this.

5.3. Charging Unifit

Unifit is powered by a built-in rechargeable Li-ion battery.

A charging adapter is included in the kit.

The battery status can be seen in the top right corner of the display.

When the remaining battery life is below 20%, the battery icon will flash. Please charge your Unifit device at this stage.

 **Warning:** Only use the charging adapter supplied. Use of other chargers could be hazardous and will negate the warranty.

To charge the battery:

Connect the charger to the unit and plug it into a mains socket. A charging screen will be displayed, showing the battery percentage level. Charge the device until the battery level is 100%.

 **Note:** The device cannot be used while being charged.

Battery life and maintenance

The battery should last at least 10 hours at 50mA, 300µs, 50Hz in TENS mode.

To fully charge the battery it will take about 2 hours.

 **Note:** The device should be charged at least once every 3 months to ensure the battery life is maintained.

 **Note:** Charging voltage: DC 5V, charging current: 90mA.

Battery warnings

- This product is equipped with a Lithium-ion battery. Failure to follow these instructions could cause the lithium-ion battery to leak acid, become hot, explode or ignite and cause injury and / or damage.
- Only use the charging adapter provided. Use of other chargers could be hazardous and will negate the warranty.
- Do not pierce, open, disassemble, or use in a humid and/or corrosive environment.
- Do not expose to temperatures over 60°C(140F).
- Do not put, store or leave near sources of heat, in direct strong sunlight, in a high temperature location, in a pressurized container or in a microwave oven.
- Do not immerse in water or get wet.
- Do not short-circuit.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with a protective earth connection.
- If battery leakage occurs and comes in contact with the skin or eyes, wash thoroughly with lots of water and immediately seek medical attention.
- Keep batteries out of the reach of children to prevent them from swallowing them by mistake. If swallowed by child, contact doctor immediately.

Battery cautions

- Do not use in a location where there is a large amount of static electricity, otherwise, the safety devices in the battery pack may be compromised.
- Disposal: Always dispose of batteries and device responsibly according to local government guidelines. Do not throw batteries onto a fire. Risk of explosion.

5.4. Electrode and lead wire care

- The electrodes are reusable but for single patient use. The electrodes should provide up to 20 uses. This may vary depending on skin condition, humidity and storage between uses.
- To maximise the life of the electrodes, always ensure the electrodes are replaced on to their plastic liner and stored and sealed in the plastic bag. When the electrodes dry, this may negatively impact the performance of the treatment.
- The adhesive in the electrode is a peelable hydrogel (water based). To restore the tackiness of the electrode, add a tiny drop of water and spread it around. Avoid adding too much, as it will make the electrode too soft. If this occurs, place the electrodes sticky side up in a refrigerator for a few hours to regain adhesiveness.
- Replace the electrodes when they lose their stickiness. Poor connection may

cause discomfort and skin irritation.

- To ensure safety and effectiveness we recommend using TensCare electrodes only.
- The electrodes provided are latex-free.
- The life of the lead wires depends greatly on use. Always handle the leads with care. We recommend replacing the lead wires regularly (about every 6 months).

5.5. Replacement parts/ accessories

The following replacement parts are compatible with the Unifit and may be ordered from TensCare at www.tenscare.co.uk or +44(0) 1372 723434.

E-CM5050	50x50mm electrodes (pack of 4)
E-CM50100	50x100mm electrodes (pack of 4)
E-CR5050	50x50mm Cool Relief infused electrodes (pack of 4)
E-CR50100	50x100mm Cool Relief infused electrodes (pack of 4)
E-TRBG	50x50mm Blue Gel hypoallergenic electrodes (pack of 4)
E-CMR60130	Breast electrodes (pack of 4)
L-CPT	Lead wire

X-BX-UR	Belt clip
X-UNICH	Charging adapter
X-MULTI-UK	UK socket
X-MULTI-EU	EU socket
X-MULTI-USA	US socket


6. Technical Specification

TENS & EMS

Amplitude	99mA zero to peak in 99 steps +-10%
Waveform	Symmetrical rectangular bi-phasic
Output	Constant current 500-1000 Ohm, Constant voltage > 1000 Ohm
Max. pulse energy	Total output limited to 25 μ C per pulse
Output plug	Fully shielded: touch proof
Channels	Dual channel
Battery	Built-in lithium-ion polymer battery Ratings: 1500mAh, 5.5Wh, 3.7V
Charging adapter	Output: DC 5 V; 1A Rated input voltage: >250 V Classification of protection against electric shock: Class II IP classification: IP21 Complies with IEC 60601-1 or IEC 60950-1
Weight	128.5g
Dimensions	125 x 62 x 28mm
Safety classification	Internal power source, designed for continuous use, no special moisture protection

Service Life	>5 years
Applied part	Skin surface electrodes
IP rating	IP22
Operating conditions	Temperature range: 5 to 40°C Humidity: 15 to 90% RH non-condensing Atmospheric pressure: 700 hPa to 1060 hPa
Storage conditions	Temperature range: -25 to 60°C Humidity: 15 to 90% RH non-condensing Atmospheric pressure: 700 hPa to 1060 hPa

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Thank you for using Unifit...

At TensCare, we strive to provide the best possible products and services. We value your feedback and continuously work to enhance our offerings. Understanding how our products benefit you is important to us. If you have any insights or suggestions, please get in touch.

Contact us

-  TensCare Ltd.
9 Blenheim Road,
Epsom, Surrey,
KT19 9BE, U.K.
-  +44 (0) 1372 723434
-  www.tenscare.co.uk
-  customer.service@tenscare.co.uk

Follow us

-  TensCare Ltd
-  @tenscareltd
-  @TensCareLtd
-  TensCare Limited

EU Declaration of Conformity

Tenscare Ltd herewith declare that the product meets the provisions of REGULATION (EU) 2017/745 which apply to it. The medical device has been assigned to class IIa according to Annex IV of the REGULATION (EU) 2017/745. The product concerned has been designed and manufactured under a quality management system according to Annex IX, of REGULATION (EU) 2017/745. This EU declaration of conformity is issued under the sole responsibility of the manufacturer. No "Common Specification" is applicable.

www.tenscare.co.uk

TensCare Ltd.

9 Blenheim Road, Epsom,
Surrey, KT19 9BE, U.K.

EU DoC

Complies with Medical Device
Regulation (EU) 2017/745.

Advena Ltd.

Tower Business Centre,
2nd Flr, Tower Street, Swatar,
BKR 4013, Malta

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