Instructions for Use

Please read carefully

Intended Use

The iTouch Go vaginal probe is a medical device designed to be used in the home healthcare environment alongside the iTouch Go Pelvic Floor Strengthener to stimulate and train the pelvic floor muscles, helping to treat symptoms of urinary incontinence. Pelvic floor stimulation also helps to relieve symptoms of chronic pelvic pain. Additionally, it may help improve sexual intimacy by toning the pelvic floor muscles.

Contraindications

Follow the Contraindications and Warnings supplied with your iTouch Go device manual.

Do not use with electrical stimulation in patients with a history of carcinoma at the site of stimulation.

Cautions

This probe is intended for single patient use only. Do not share your probe with anyone else. Improper treatment or cross-infection may occur.

It is important that the probe is cleaned before and after use. Ineffective cleaning may lead to irritation or infection.

Never insert or remove the probe unless the control unit is

powered OFF as insertion or removal when stimulation is active may cause discomfort or tissue irritation.

If tissue irritation occurs, discontinue treatment immediately.

Ask your healthcare professional for advice before continuing treatment to prevent injury.

Do not use a silicone-based lubricant on the metal plates of the probe as it may decrease the effectiveness of the electrical stimulation.

Electrode surface area is 6.08cm². At high output settings on some devices, this could result in a current density greater than 2mA RMS/cm². High current density may be associated with localised tissue damage or minor non-thermal burns. Increase intensity carefully and reduce if you experience sharp pain during use or discharge or bleeding after use. Do not use if the pack is opened or damaged.

How to Use

Connecting the replacement probe

1. Ensure the iTouch Go control hub is turned off and remove the old probe.
2. Place the new probe in the control hub.
3. Connect the control hub to a power source via the USB charging cable.
4. Wait until the probe starts charging. This will be indicated by the probe battery icon and green light flashing.
5. The new probe is now connected to the control hub and ready for use.

Warning: Only connect one probe to the control hub at any one time.

Using the probe

1. Clean the probe (see cleaning instructions).
2. Before using the probe, you will need to visit the toilet.
3. Lubricate the probe with a water-based lubricant, such as TensCare Go Gel.
4. Choose a comfortable position, such as lying down on your bed on your side with your knees raised.
5. Insert the probe into your vagina, until only the plastic rim at the base of the probe is visible.
6. Follow guidance on completing a session in the iTouch Go device instructions for use.

Caution: Do not use a silicone-based lubricant on the stimulation contacts as it may decrease the effectiveness of the muscle stimulation.

Warning: Ensure the iTouch Go is switched off before inserting the probe.

After a session

1. Check the iTouch Go control hub is off. If it is not, hold down the power button on the control hub to switch off.

itouch

Vaginal Probe go

2. Remove the probe from the vagina by holding the probe tail and gently pulling outwards.
3. Clean and dry the probe before storing it in the iTouch Go control hub.

Charging

The iTouch Go control hub acts as a charging hub for the wireless probe. The hub features a Type C charging socket found on the side of the control hub.

1. Insert the probe into the control hub.
2. Connect the control hub to a power source using the charging cable provided in the iTouch Go Kit.
3. The iTouch Go probe battery level is displayed on the top right of the iTouch Go control hub screen. The LED indicator below the probe flashes green when the probe is charging and remains lit when fully charged.

Note: The probe should be charged at least once every 3 months to ensure battery life is maintained.

Battery warnings

This product is equipped with a Lithium-ion battery. Failure to follow these instructions could cause the lithium-ion battery to leak acid, become hot, explode, or ignite and cause injury and/or damage:

Do not pierce, open, disassemble, or use in a humid and/or corrosive environment.

Do not expose to temperatures over 60°C(140°F).

Do not put, store, or leave near sources of heat, in direct strong sunlight, in a high temperature location, in a pressurized container or in a microwave oven.

Do not immerse in water.

If battery leakage occurs and comes in contact with the skin or eyes, wash thoroughly with lots of water and immediately seek medical attention.

Do not use in a location where there is a large amount of static electricity, otherwise, the safety devices in the battery pack may be compromised.

Disposal: Always dispose of batteries and device responsibly according to local government guidelines. Do not throw batteries onto a fire. Risk of explosion.

Cleaning

It is important that the probe is cleaned before and after each use. Clean with either an alcohol-free antibacterial wipe such as TensCare Wipes (see X-WIPES) or by wiping with warm soapy water. Rinse under running water and dry thoroughly and return the unit to the storage pouch.

Warning: Do not immerse the probe in a liquid.

Troubleshooting

If you are using the probe, you are experiencing no sensation and the display on the control hub is flashing, the probe may be faulty. Test the probe by holding the probe in your hand.

1. Dampen your hand with water and a little table salt.
2. Squeeze the probe firmly and make sure your skin is covering the metal parts of the probe and carefully increase strength until you can feel something.
3. Most people will feel the stimulation at around 25 mA.
4. If the display flashes and the unit will not allow you to pass 20 mA, the probe is faulty and needs replacing.

To report unexpected operation or events contact Tenscare.

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