

## Quickstart guide

### Preparation before treatment:

1) Wash your face (or area being treated) with clear water and dry it.



2) Attach the light module with the control unit. Ensure that the golden connectors from the light module are in contact with the golden pins from the control unit.



3) Insert the light gathering cover (if not already on) to ensure the light is 2 cm away from the skin. When using the light gathering cover, you can simply use the unit with the light gathering cover being in contact with the skin.



### Control buttons:

Button	Description
	Power switch
	Work mode
	Start/stop
	Start/stop pulse mode

### Treatment session:

- WARNING:** Do not direct the lights of the unit towards your eyes. You can simply direct the light towards your hand while you are setting up the unit.
- Press and hold down the ON/OFF button for 3 seconds, , at the top back of the unit. The unit will automatically enter the last mode used.
- Press the work mode button, , to choose a different work mode. You can choose 1 (1 min), 2 (2 min), 3 (3 min), 5 (5 min) or L (continuous).
- Press the pulse mode button, , if you want to activate or deactivate the pulse mode. Using the pulse mode may help cell regeneration (anti-ageing). Using the normal mode may help speed up healing (acne treatment).
- Cover both eyes with the special goggles provided to avoid light source damaging them.
- Now place the unit, with the light gathering cover in contact with the skin, onto the area you want to treat and press the start/stop button, .
- If you selected a time, the unit will stop emitting light and beep when the treatment time is over. If you have selected L (Continuous), press the start/stop button to stop the emission of light.

- Then, simply reposition the unit to the next area to treat and press the start/stop button again.
- Once the desired areas have been treated, you can switch off the unit by pressing and holding down the ON/OFF button for 3 seconds, .

### Treatment time and frequency:

Module	Light	Treatment time	Indication
For wrinkles	Red and infrared* light	2 to 3 minutes daily use. Once good results have been achieved, use 2 to 3 times per week.	Diminish fine lines, firm skin, increase collagen, increase the skin's moisture. Fully absorbed by fibre cells to promote cell growth; stimulate cells to produce collagen. Thickening and restructuring the dermis structure, smooth and increase the elasticity of the skin.
For acne improvement	Blue and red lights	5 minutes daily use. Once good results have been achieved, use 2 to 3 times per week.	Treat acne and repair damaged skin. Activate the production of the protoporphyrin within acne, leads to single oxygen release, kill Propionibacterium within acne, rapidly remove pustule, improve the skin surface.

\*Please note that the infrared lights are not visible.

### Troubleshooting:

Problem	Possible causes	Solution
LED is not working properly	Low power	Charge the unit using the USB charging lead
	LED light module is not connected properly	Reposition the module
	Battery installation is not correct	Reposition the battery and check the polarities
	Broken battery	Contact customer service
Cannot get it to start	Battery installation is not correct	Reposition the battery and check the polarities
	The control unit is damaged/faulty	Contact customer service
	Broken battery	Contact customer service
Cannot charge	Battery installation is not correct	Reposition the battery and check the polarities
	The USB charging lead is damaged/faulty	Contact customer service
	Broken battery	Contact customer service

## Quickstart guide

### Preparation before treatment:

1) Wash your face (or area being treated) with clear water and dry it.



2) Attach the light module with the control unit. Ensure that the golden connectors from the light module are in contact with the golden pins from the control unit.



3) Insert the light gathering cover (if not already on) to ensure the light is 2 cm away from the skin. When using the light gathering cover, you can simply use the unit with the light gathering cover being in contact with the skin.



### Control buttons:

Button	Description
	Power switch
	Work mode
	Start/stop
	Start/stop pulse mode

### Treatment session:

- WARNING:** Do not direct the lights of the unit towards your eyes. You can simply direct the light towards your hand while you are setting up the unit.
- Press and hold down the ON/OFF button for 3 seconds, , at the top back of the unit. The unit will automatically enter the last mode used.
- Press the work mode button, , to choose a different work mode. You can choose 1 (1 min), 2 (2 min), 3 (3 min), 5 (5 min) or L (continuous).
- Press the pulse mode button, , if you want to activate or deactivate the pulse mode. Using the pulse mode may help cell regeneration (anti-ageing). Using the normal mode may help speed up healing (acne treatment).
- Cover both eyes with the special goggles provided to avoid light source damaging them.
- Now place the unit, with the light gathering cover in contact with the skin, onto the area you want to treat and press the start/stop button, .
- If you selected a time, the unit will stop emitting light and beep when the treatment time is over. If you have selected L (Continuous), press the start/stop button to stop the emission of light.

- Then, simply reposition the unit to the next area to treat and press the start/stop button again.
- Once the desired areas have been treated, you can switch off the unit by pressing and holding down the ON/OFF button for 3 seconds, .

### Treatment time and frequency:

Module	Light	Treatment time	Indication
For wrinkles	Red and infrared* light	2 to 3 minutes daily use. Once good results have been achieved, use 2 to 3 times per week.	Diminish fine lines, firm skin, increase collagen, increase the skin's moisture. Fully absorbed by fibre cells to promote cell growth; stimulate cells to produce collagen. Thickening and restructuring the dermis structure, smooth and increase the elasticity of the skin.
For acne improvement	Blue and red lights	5 minutes daily use. Once good results have been achieved, use 2 to 3 times per week.	Treat acne and repair damaged skin. Activate the production of the protoporphyrin within acne, leads to single oxygen release, kill Propionibacterium within acne, rapidly remove pustule, improve the skin surface.

\*Please note that the infrared lights are not visible.

### Troubleshooting:

Problem	Possible causes	Solution
LED is not working properly	Low power	Charge the unit using the USB charging lead
	LED light module is not connected properly	Reposition the module
	Battery installation is not correct	Reposition the battery and check the polarities
	Broken battery	Contact customer service
Cannot get it to start	Battery installation is not correct	Reposition the battery and check the polarities
	The control unit is damaged/faulty	Contact customer service
	Broken battery	Contact customer service
Cannot charge	Battery installation is not correct	Reposition the battery and check the polarities
	The USB charging lead is damaged/faulty	Contact customer service
	Broken battery	Contact customer service